

## MultiMini-BH Series Multi Mini Block Heaters

**MultiMini-BH Series** are a dry bath incubator with 3 different block heating features at the same time. The BH3 is the incubator that has 3 individually controlled dry mini-baths and The BH3C is incubator that has for both heating and cooling feature. Compatible with many aluminum blocks covering various tube volumes. This series is suitable for laboratories with heavy workloads.

### Features;

1. 3 independent-controlled block heater systems that save working space, time and budget.
2. Easy to Install and use LCD Screen.
3. Simultaneous preview of set time and temperature values and real time and temperature values.
4. Automatic stop and alarm after program completion.
5. Controllable heating and cooling system from 0 ° C to 99 ° C. (For BH3C/Block-A)

### Technical Specifications;

Model	MultiMini-BH3C			MultiMini-BH3
Working Area	Block-A	Block-B	Block-C	Block-A/B/C
Temperature Control Range	0°C - 99.9°C	+5 - 99.9°C	+5 - 99.9°C	+5 - 99.9°C
Heating Time (25 to 99.9°C)	≤20min	≤20min	≤20min	≤20min
Cooling Time (25 to 0°C)	≤25min (*)	N/A	N/A	N/A
Temperature Control Accuracy	≤0.5°C			≤0.5°C
Temp. Control Accuracy at 37°C	0.2°C			0.2°C
Display Accuracy	0.1°C			0.1°C
Timing Range	1 - 999h			1 - 999h
Block Quantity	3			3
Power	DC12V,120W			DC12V, 120W
Dimension(mm)	360x180x150			360x180x150
Net Weight	3kg			3kg



### Block Types;

Type	Description
A-MiniT	40x0.2ml
B-MiniT	24x0.5ml
C-MiniT	15x1.5ml
D-MiniT	15x2.0ml
E-MiniT	8--12.5x12.5x32mm
F-MiniT	4x15ml
G-MiniT	2x50ml

