

INOCI Series High Capacity Cooling Incubators



INOCI Series High Capacity Refrigerating Incubators are made of high quality stainless steel plate and it has frozen cold storage function.

With its easy to clean internal structure, it reduces the need for work force and provides easy operation with simple interface.

R134a Refrigerant gas offers healthy and reliable operation.

CE



Technical Specifications;

Model	INOCI-C1020	INOCI-SC1020
Capacity		1020L
Temperature Range	0 - 10°C	-10 - 60°C
Refrigerant		R134a
Cooling Mode		Fan cooling or Static
Chamber Material		Stainless Steel
Power Supply		220V 50Hz
Power Consumption		1500W
Interior Dimension		1100x630x1470mm
Exterior Dimension		1220x780x2050mm
Weight		161 Kg