

## INOWB Series Multi Water Bath Systems

INOWB Series Water Baths are apply to distillate, drying, concentrate and thermostatic heat chemical drugs, biological products etc in those laboratory department such as medical treatment units, scientific researches, chemical industry, pharmacy, oil, printing, painting and etc.

Thanks to the inner chamber and lid made of stainless steel, it provides a long and healthy use. Intelligent structure design with strong corrosion resistance, strong appearance and easy to use.

With the over-temperature alarm system, it ensures that the samples are safely operated.

The microcomputer makes it easy to use, providing temperature control, high precision temperature control and stable performance.



CE



### Features;

1. Extreme temperature alarm system.
2. Microcomputer PID temperature control.
3. High-powered CPU chip.
4. High precision and precision platinum resistor PT100 sensor.
5. Existing timing function.

### Technical Specifications;

Model	INOWB-1	INOWB-2	INOWB-4	INOWB-6	INOWB-8
Number of Chamber	1	2	4	6	8
Capacity (L)	2.5	4.9	9.9	14.8	19.8
Temperature Range	RT +5°C - 99°C				
Timer Range	1 - 999min				
Temperature Fluctuation	±0.2°C				
Temperature Tracking Alarm	±2°C				
Power Supply	220V ±10% 50/60Hz				
Power Consumption (W)	300	500	1000	1500	2000
Dimensions (LxWxH)(mm)	350x260x220	570x240x270	670x450x260	800x450x270	800x450x270
Package Size (LxWxH)(mm)	450x360x320	670x340x370	770x550x360	900x550x370	900x550x370
Weight (kg)	7.5	12.5	16	20.5	35