

INOOB Series Electric Thermostatic Oil Bath

INOOB Series Oil Baths, Anticorrosive stainless steel and microcomputer operating system provide the best safety and balance.

Features;

1. New design and ergonomic design.
2. Stainless steel hopper and lid and oven polish applied surface.
3. The microprocessor PID control with timing function, high-power CPU chip and high-precision platinum resistance PT100 sensor, makes the temperature more accurate and easier to operate.
4. The microcomputer controls the speed of the mixture to ensure that the rotation speed is constant under the constantly changing mixing viscosity condition.
- * Only Magnetic Stirrer Bath has this function.
5. Mixer function ensures more precise and uniform temperature control.



Technical Specifications;

Model	INOOB-5	INOOB-13	INOOB-28
Capacity (L)	5	13	28
Temperature Range	RT +5°C - 200°C		
Timer Range	1 - 999min		
Temperature Stability	±0.5°C		
Tracking Alarm	±2°C		
Power Supply	AC220V 50/60Hz		
Power Consumption (W)	500	800	1800
Inner Dimension(LxWxH)(mm)	300x150x150	295x235x200	495x295x200
Dimensions (LxWxH)(mm)	510x250x270	940x480x270	940x480x270
Package Size (LxWxH)(mm)	600x350x370	1040x580x370	1040x580x370
Weight (kg)	9.3	15.8	23