

INOMINO Series Mini Centrifuges

INOMINO Series Mini Centrifuges are supplied 6 x 1.5ml rotor for micro tubes.

The 6 x 1.5ml rotor is designed to centrifuge up to six individual 1.5ml plastic micro centrifuge tubes.

It will also accept 0.5ml and 0.2ml tubes with the adapters supplied with the unit.

Adapters are available separately for 0.2ml tubes and 0.5ml tubes.

The rotor is designed for applications requiring relatively low g-forces, such as micro filtration, cell separation and quick spin downs from the walls of tubes.

When the appliance is connected to electricity, centrifugation starts when the lid is closed. The sensor that detects the lid condition also automatically stops the operation when the lid is opened.

Features;

The device is equipped with a detailed and colorful cover for monitoring at runtime.

The rotor it contains is made of Carbon Fiber which is a high quality material.

The INOMINO series occupies a small area and operates at low noise.

It has all the features that require it to be safe and durable.



Technical Specifications;

Model	INOMINO-6K	INOMINO-7K	INOMINO-10K
Max. RPM	6000	7000	10000
Relative Centrifugal Force	2000	2300	5000
Rotor Types	6 x 1.5/2.0ml Angle Rotor 2 x 8 x 0.2ml Strip Rotor		
Power Supply	30/40W, AC100V-240V, 50/60Hz, 1A		
Dimension	175x150x120mm		
Packing Dimension	195x260x160mm		
Weight	0,7 Kg		
Packing Weight	2,5 Kg		

