


BRC Series Bench-Top Touch Panel Refrigerated Centrifuges Rotor Types



MODEL		BRC-MICRO		BRC-5180		BRC-5300		BRC-5800	
	Capacity	RPM	RCF(xg)	RPM	RCF(xg)	RPM	RCF(xg)	RPM	RCF(xg)
Angel Rotor	40x0.2ml	16000	19041						
	24x0.5ml	16000	18481						
	12x1.5/2ml	16000	17948						
	10x5ml	16000	17882						
	20x1.5/2ml	14000	15585						
	24x1.5/2ml	14000	17956						
	6x10/15ml	13000	14196						
	40x0.5ml	13000	17218						
	30x0.5ml	13000	13906						
	4x8 PCR	14000	12720						
	24 Hematocrit Rotor	12000	15805						
	16x10/15ml			4000	2820				
	24x10/15ml			4000	2820				
	6x15ml			5000	2540			5000	2540
	12x15ml			5000	3080	5000	3080	5000	3080
	24x15ml			5000	3500	5000	3500	5000	3500
	30x15ml							5000	3830
	4x50ml			5000	2520	5000	2520	5000	2520
	6x50ml			5000	2850	5000	2850	5000	2850
	12x50ml					5000	3860	5000	3860
	24x50ml							5000	2970
Swing Rotor	4x100ml			5000	2630	5000	2630	5000	2630
	6x100ml			5000	3130	5000	3130	5000	3130
	12x100ml					4000	2970	4000	2970
	4x180ml			4000	2580				
	2x2x96			4000	2020				
	4x50ml			5000	4420				
	4x300ml(Square Cup)					4000	2920		
	4x300ml(Round Cup)					4000	2920		
	4x500ml(Square Cup)							4000	3310
	4x500ml(Round Cup)							4000	3380
	4x800ml(Round Cup)							4000	3450
	4x24x5ml					4000	2840		
	4x24x7ml					4000	3140		
	4x12x10/15ml					4000	3140		
	4x31(Khan Tube)					4000	2840		
	4x50ml							5000	4730
	8x50ml							4000	3200
	4x100ml							5000	4730
	8x100ml							4000	3020
	16x15ml							4000	3200
	24x15ml							4000	3200
	32x15ml							4000	3200
	48x10/15ml							4000	2880
	48x7ml							4000	2760
	48x5ml							4000	2480
	72x10/15ml							4000	3140
	88x7ml							4000	3140
	88x5ml							4000	2840
	16x10/7/5ml							4000	2790
	24x10/7/5ml							4000	2790
	32x10/7/5ml							4000	2790
Microplate Rotor	2x4x96 Well					4000	2490	4000	2490
	2x3x48 Well					4000	2300	4000	2300
	4x4x96 Well							4000	2840